


## THE MOST AGILE SURFACE TEMPERATURE SENSOR

eTemp Performance is a unique solution that continuously collects, records and transmits **surface temperature** thanks to **Bluetooth Low Energy** communication.

High accuracy:  $\pm 0.1^{\circ}\text{C}$



Waterproof and super resistant  
to use in **harsh conditions**



- **Rechargeable,**
- Real time and a *posteriori* data collection,
- **Easy to use.**

- Follow **multiple sensors** at a time,
- Get access to raw data in **.CSV** or **PDF**.

### Technical Features:

- **Temperature range:**  $-30^{\circ}\text{C}$  to  $50^{\circ}\text{C}$
- **Temperature accuracy:**  $\pm 0.1^{\circ}\text{C}$
- **Temperature resolution:**  $0.01^{\circ}\text{C}$
- **Dimensions:** 12.4 x 9.5 x 5.8 mm
- **Weight:**  $\approx 1.3$  g
- **Memory capacity:** 3600 temperature values can be stored into the sensor's internal memory
- **1 kit:** 4 sensors and 1 charger box



BodyCAP  
Your e-health partner